

Irrigation
solutions

Calix



Calix

SODIUM HYPOCHLORITE IRRIGATION SOLUTIONS

- Ideal for debridement and irrigation of root canals during and after instrumentation
- Gluten-Free
- Made in USA

3% SODIUM HYPOCHLORITE SOLUTION

- Item #28-00001
3% Sodium Hypochlorite Solution - 125ml bottle
- Item #28-00002
3% Sodium Hypochlorite Solution - 500ml bottle



6% SODIUM HYPOCHLORITE SOLUTION

- Item #28-00006
6% Sodium Hypochlorite Solution - 125ml bottle
- Item #28-00007
6% Sodium Hypochlorite solution - 500ml bottle



Calix^X

6% ENHANCED SODIUM HYPOCHLORITE SOLUTION

- Item #28-00011
6% Enhanced Sodium Hypochlorite Solution
125ml bottle
- Item #28-00012
6% Enhanced Sodium Hypochlorite Solution
500ml bottle



- Calix-X is an enhanced 6% sodium hypochlorite solution designed for irrigation, debridement and cleaning of root canals during and after instrumentation.
- Contains wetting agents and surface modifiers that enables it to penetrate into lateral canals and isthmuses by lowering surface tension.
- Gluten-Free
- Reduces post-treatment complications. Improves contact and dissolves tissue.
- Increases lubricity and decreases surface tension.
- Made in USA

Calix^E

17% EDTA ROOT CANAL SOLUTION

- Item #28-00020
17% EDTA Solution – 125ml bottle
- Item #28-00021
17% EDTA Solution – 500ml bottle



- EDTA – Ehtylenediaminetetraacetic acid in an aqueous solution.
- Calix-X is a root canal chelating agent that conditions and cleans through a chelation process. It is the irrigant of choice for smear layer and dentin mud and can be used as a final irrigant prior to obturation.
- It is perfect as a lubricant when instrumenting tight canals because it contains a surfactant which enables the EDTA to wet the canal easily. It prepares dentin walls for better adhesion of filling materials.
- It opens dentin tubules and enhances cleansing.
- It is buffered to an optimum neutral PH. For best outcomes, the PH is more important than the concentration. As PH increases, the availability of calcium ions from hydroxyapatite for chelation decreases.
- Gluten Free
- Made in USA

Calix^{CHG}

2% CHLORHEXIDINE GLUCONATE SOLUTION

- Item #28-00025
2% CHX Solution – 125ml bottle
- Item #28-00026
2% CHX Solution – 500ml bottle



- Calix-CHG contains 2% Chlorhexidine Gluconate plus wetting agents and surface modifiers to improve the depth of penetration.
- It kills bacteria faster than standard 2% Chlorhexidine solutions.
- Very effective in killing the E faecalis bacteria. It has a wide antimicrobial spectrum and is effective against both Gram-negative and Gram-positive bacteria as well as yeasts.
- Gluten-Free
- Made in USA

Calix^R

95% ETHYL ALCOHOL RINSE

- Item #28-00031
95% Ethyl Alcohol Solution
125ml bottle



- Ideal rinse for penetrating lateral canals and isthmuses to improve drying and disinfection. It has a low surface tension.
- Ideal for use after the smear layer removal with Calix-E.
- It clears the canal of all sodium hypochlorite before the final irrigation with Calix-CHG.
- Gluten-Free
- Made in USA

Distributed by: