

jade

ETCHANT GELS



Jade is a **37% thixotropic phosphoric, non-slump acid gel** for enamel and dentin. Jade is for those clinicians who prefer a creamy gel that flows into dentin and enamel, but does not run.

- Jade stays where is placed and **does not run** onto soft tissue.
- Jade **washes off completely** and easily and **does not stain**.
- Jade has a dark blue color that allows **high visibility** for **precise selective placement**. It is also available in a lighter green.
- Jade creates micropores in the enamel surface and removes the smear layer to unblock dentinal tubules.
- Gluten-Free
- Made in USA

BLUE

- Item #20-00001
Uni-Syringe Econo Kit
1 x 12g syringe
5 Pre-Bent 23ga Blue Applicator Tips
- Item #20-00002
Uni-Syringe Kit
1 x 12g syringe
20 Pre-Bent 23ga Blue Applicator Tips
- Item #20-00003
Mini Kit
4 x 1.2ml syringes
20 Pre-Bent 23ga Blue Applicator Tips



- Item #20-00004
Mini Bulk Kit
25 x 1.2ml syringes
50 Pre-Bent 23ga Blue Applicator Tips



- Item #20-00005
Bulk Kit
1 x 50ml syringe, 4 x 1.2ml empty syringes
50 Pre-Bent 23ga Blue Applicator Tips

- Item #20-00016
1 x 50ml
refill syringe



GREEN

- Item #20-00006
Uni-Syringe Econo Kit
1 x 12g syringe
5 Pre-Bent 23ga Blue Applicator Tips
- Item #20-00007
Uni-Syringe Kit
1 x 12g syringe
20 Pre-Bent 23ga Blue Applicator Tips
- Item #20-00008
Mini Kit
4 x 1.2ml syringes
20 Pre-Bent 23ga Blue Applicator Tips



- Item #20-00009
Mini Bulk Kit
25 x 1.2ml syringes
50 Pre-Bent 23ga Blue Applicator Tips



- Item #20-00010
Bulk Kit
1 x 50ml syringe, 4 x 1.2ml empty syringes
50 Pre-Bent 23ga Blue Applicator Tips

- Item #20-00091
1 x 50ml
refill syringe



APPLICATORS

- Item #20-00011 - 25 Jade Pre-Bent 23ga Blue Applicator Tips
- Item #20-00012 - 100 Jade Pre-Bent 23ga Blue Applicator Tips
- Item #20-00013 - 25 Jade Pre-Bent 25ga Blue Applicator Tips
- Item #20-00014 - 100 Jade Pre-Bent 25ga Blue Applicator Tips

Distributed by: